

# ITPK1 Antibody



PACO10017

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ITPK1

**Uniprot**

Q13572

**Synonyms:**

inositol 1,3,4-triphosphate 5/6 kinase; ITPK1; ITRPK1 ;

**Immunogen:**

Human ITPK1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A