

ITPKC Antibody



PACO10020

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ITPKC

Uniprot

Q96DU7

Synonyms:

inositol 1,4,5-trisphosphate 3-kinase C; ITPKC; IP3KC ;

Immunogen:

Human ITPKC.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A