

IVNS1ABP Antibody



PACO10022

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

IVNS1ABP

Uniprot

Q9Y6Y0

Synonyms:

influenza virus NS1A binding protein; IVNS1ABP; DKFZp686K06216; FLARA3; FLJ10069; FLJ10411; FLJ10962; FLJ35593; FLJ36593; HSPC068; KIAA0850; ND1; NS-1; NS1-BP; NS1BP ;

Immunogen:

Human IVNS1ABP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A