

JMJD6 Antibody



PACO10028

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

JMJD6

Uniprot

Q6NYC1

Synonyms:

jumonji domain containing 6; JMJD6; KIAA0585; PSR; PTDSR; PTDSR1 ;

Immunogen:

Human JMJD6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A