

JSRP1 Antibody



PACO10030

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

JSRP1

Uniprot

Q96MG2

Synonyms:

junctional sarcoplasmic reticulum protein 1; JSRP1; FLJ32416; JP-45 ;

Immunogen:

Human JSRP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A