

PACO10037

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KARS

Uniprot

Q15046

Synonyms:

lysyl-tRNA synthetase; KARS; KARS2; KIAA0070; KRS ;

Immunogen:

Human KARS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A