

KAT5 Antibody



PACO10040

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KAT5

Uniprot

Q92993

Synonyms:

K(lysine) acetyltransferase 5; KAT5; ESA1; HTATIP; HTATIP1; PLIP; TIP; TIP60; cPLA2 ;

Immunogen:

Human KAT5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A