

KBTBD2 Antibody



PACO10043

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KBTBD2

Uniprot

Q8IY47

Synonyms:

kelch repeat and BTB (POZ) domain containing 2; KBTBD2; BKLHD1 ;

Immunogen:

Human KBTBD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A