

KCNA2 Antibody



PACO10044

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCNA2

Uniprot

P16389

Synonyms:

potassium voltage-gated channel, shaker-related subfamily, member 2; KCNA2; HBK5; HK4; HUKIV; KV1.2; MGC50217; MK2; NGK1; RBK2 ;

Immunogen:

Human KCNA2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A