

KCND2 Antibody



PACO10048

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCND2

Uniprot

Q9NZV8

Synonyms:

potassium voltage-gated channel, Shal-related subfamily, member 2; KCND2; KIAA1044; KV4.2; MGC119702; MGC119703; RK5 ;

Immunogen:

Human KCND2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A