

KCNH7 Antibody



PACO10052

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KCNH7

Uniprot

Q9NS40

Synonyms:

potassium voltage-gated channel, subfamily H (eag-related), member 7; KCNH7; ERG3; HERG3; Kv11.3; MGC45986 ;

Immunogen:

Human KCNH7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A