

# KCNJ10 Antibody



PACO10057

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KCNJ10

**Uniprot**

P78508

**Synonyms:**

potassium inwardly-rectifying channel, subfamily J, member 10; KCNJ10; BIRK-10; KCNJ13-PEN; KIR1.2; KIR4.1 ;

**Immunogen:**

Human KCNJ10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A