

KCNJ15 Antibody



PACO10061

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KCNJ15

Uniprot

Q99712

Synonyms:

potassium inwardly-rectifying channel, subfamily J, member 15; KCNJ15; IRKK; KIR1.3; KIR4.2; MGC13584 ;

Immunogen:

Human KCNJ15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A