

KCNJ6 Antibody



PACO10064

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCNJ6

Uniprot

P48051

Synonyms:

potassium inwardly-rectifying channel, subfamily J, member 6; KCNJ6; BIR1; GIRK2; KATP2; KCNJ7; KIR3.2; MGC126596; hiGIRK2 ;

Immunogen:

Human KCNJ6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A