

KCNK5 Antibody



PACO10066

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KCNK5

Uniprot

O95279

Synonyms:

potassium channel, subfamily K, member 5; KCNK5; FLJ11035; K2p5.1; TASK-2; TASK2 ;

Immunogen:

Human KCNK5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A