

KCTD15 Antibody



PACO10071

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

KCTD15

Uniprot

Q96SI1

Synonyms:

potassium channel tetramerisation domain containing 15; KCTD15; MGC25497; MGC2628 ;

Immunogen:

Human KCTD15.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A