

KDM3B Antibody



PACO10073

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KDM3B

Uniprot

Q7LBC6

Synonyms:

lysine (K)-specific demethylase 3B; KDM3B;5qNCA; C5orf7; JMJD1B; KIAA1082 ;

Immunogen:

Human KDM3B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A