

KDM5C Antibody



PACO10075

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

KDM5C

Uniprot

P41229

Synonyms:

lysine (K)-specific demethylase 5C; KDM5C; DXS1272E; JARID1C; MRXJ; MRXSJ; SMCX; XE169 ;

Immunogen:

Human KDM5C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A