

KIAA0391 Antibody



PACO10083

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

KIAA0391

Uniprot

O15091

Synonyms:

KIAA0391; KIAA0391; MRPP3 ;

Immunogen:

Human KIAA0391.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A