

KIF26B Antibody



PACO10100

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KIF26B

Uniprot

Q2KJY2

Synonyms:

kinesin family member 26B; KIF26B; FLJ10157; MGC35030 ;

Immunogen:

Human KIF26B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A