

# KIR2DS4 Antibody



PACO10113

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KIR2DS4

**Uniprot**

P43632

**Synonyms:**

killer cell immunoglobulin-like receptor, two domains, short cytoplasmic tail, 4;  
KIR2DS4; CD158I; KIR1D; KIR412; KKA3; MGC120019; MGC125315; MGC125317; NKAT8  
;

**Immunogen:**

Human KIR2DS4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A