

KIRREL3 Antibody



PACO10115

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KIRREL3

Uniprot

Q8IZU9

Synonyms:

kin of IRRE like 3 (Drosophila); KIRREL3; KIAA1867; KIRRE; MGC126824; MGC126850; MRD4; NEPH2; PRO4502 ;

Immunogen:

Human KIRREL3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A