

KLC2 Antibody



PACO10119

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KLC2

Uniprot

Q9H0B6

Synonyms:

kinesin light chain 2; KLC2; FLJ12387 ;

Immunogen:

Human KLC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A