

KLC3 Antibody



PACO10120

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

KLC3

Uniprot

Q6P597

Synonyms:

kinesin light chain 3; KLC3; KLC2; KLC2L; KLCt; KNS2B ;

Immunogen:

Human KLC3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A