

KLF15 Antibody



PACO10123

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

KLF15

Uniprot

Q9UIH9

Synonyms:

Kruppel-like factor 15; KLF15; DKFZp779M1320; KKLF ;

Immunogen:

Human KLF15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A