

# KLHL29 Antibody



PACO10137

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

KLHL29

**Uniprot**

Q96CT2

**Synonyms:**

kelch-like 29 (Drosophila); KLHL29; KBTBD9; KIAA1921 ;

**Immunogen:**

Human KLHL29.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A