

KLHL32 Antibody



PACO10139

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

Protein Background:

Gene ID:

KLHL32

Uniprot

Q96NJ5

Synonyms:

470aa,53kd;

Immunogen:

Human KLHL32.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A