

KNG1 Antibody



PACO10150

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KNG1

Uniprot

P01042

Synonyms:

kininogen 1; KNG1; BDK; KNG ;

Immunogen:

Human KNG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A