

KPNA2 Antibody



PACO10152

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KPNA2

Uniprot

P52292

Synonyms:

karyopherin alpha 2 (RAG cohort 1, importin alpha 1); KPNA2; IPOA1; QIP2; RCH1; SRP1alpha ;

Immunogen:

Human KPNA2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A