

KPNB1 Antibody



PACO10156

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

KPNB1

Uniprot

Q14974

Synonyms:

karyopherin (importin) beta 1; KPNB1; IMB1; IPO1; IPOB; Impnb; MGC2155; MGC2156; MGC2157; NTF97 ;

Immunogen:

Human KPNB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A