

KRT19 Antibody



PACO10166

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KRT19

Uniprot

P08727

Synonyms:

keratin 19; KRT19; CK19; K19; K1CS; MGC15366 ; Cytokeratin 19;

Immunogen:

Human KRT19.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A