

KRT20 Antibody



PACO10167

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

KRT20

Uniprot

P35900

Synonyms:

keratin 20; KRT20; CD20; CK20; K20; KRT21; MGC35423 ;

Immunogen:

Human KRT20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A