

LACTB2 Antibody



PACO10184

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

LACTB2

Uniprot

Q53H82

Synonyms:

lactamase, beta 2; LACTB2; CGI-83 ;

Immunogen:

Human LACTB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A