

LAS1L Antibody



PACO10199

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LAS1L

Uniprot

Q9Y4W2

Synonyms:

LAS1-like (*S. cerevisiae*); LAS1L; FLJ12525; dJ475B7.2 ;

Immunogen:

Human LAS1L.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A