

PACO10207

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LBR

Uniprot

Q14739

Synonyms:

lamin B receptor; LBR; DHCR14B; FLJ43126; LMN2R; MGC9041; PHA ;

Immunogen:

Human LBR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A