

LDHAL6A Antibody



PACO10218

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

LDHAL6A

Uniprot

Q6ZMR3

Synonyms:

lactate dehydrogenase A-like 6A; LDHAL6A; MGC23940 ;

Immunogen:

Human LDHAL6A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A