

LECT2 Antibody



PACO10225

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LECT2

Uniprot

O14960

Synonyms:

leukocyte cell-derived chemotaxin 2; LECT2; MGC126628; chm-II; chm2 ;

Immunogen:

Human LECT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A