

LEO1 Antibody



PACO10231

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LEO1

Uniprot

Q8WVCO

Synonyms:

Leo1, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae); LEO1; RDL ;

Immunogen:

Human LEO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A