

LEPREL1 Antibody



PACO10233

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

P3H2

Uniprot

Q8IVL5

Synonyms:

leprecan-like 1; LEPREL1; FLJ10718; MLAT4; P3H2 ;

Immunogen:

Human LEPREL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A