

LGALS2 Antibody



PACO10236

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

LGALS2

Uniprot

P05162

Synonyms:

lectin, galactoside-binding, soluble, 2; LGALS2; HL14; MGC75071 ;

Immunogen:

Human LGALS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A