

LGI2 Antibody



PACO10241

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LGI2

Uniprot

Q8N0V4

Synonyms:

leucine-rich repeat LGI family, member 2; LGI2; FLJ10675; KIAA1916; LGIL2; MGC126808; MGC126810 ;

Immunogen:

Human LGI2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A