

LHCGR Antibody



PACO10245

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LHCGR

Uniprot

P22888

Synonyms:

lutinizing hormone/choriagonadotropin receptor; LHCGR; FLJ41504; LCGR; LGR2; LH/CG-R; LH/CGR; LHR; LHRHR; LSH-R;

Immunogen:

Human LHCGR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A