

LHX5 Antibody



PACO10246

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LHX5

Uniprot

Q9H2C1

Synonyms:

LIM homeobox 5; LHX5; MGC129689 ;

Immunogen:

Human LHX5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A