

LIMS1 Antibody



PACO10258

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

LIMS1

Uniprot

P48059

Synonyms:

LIM and senescent cell antigen-like domains 1; LIMS1; PINCH; PINCH1 ;

Immunogen:

Human LIMS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A