

LIN7B Antibody



PACO10260

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LIN7B

Uniprot

Q9HAP6

Synonyms:

lin-7 homolog B (C. elegans); LIN7B; LIN-7B; MALS-2; MALS2; VELI2 ;

Immunogen:

Human LIN7B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A