

LINGO1 Antibody



PACO10263

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LINGO1

Uniprot

Q96FE5

Synonyms:

leucine rich repeat and Ig domain containing 1; LINGO1; FLJ14594; LERN1; LRRN6A; MGC17422; UNQ201 ;

Immunogen:

Human LINGO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A