

# LMBR1L Antibody



PACO10268

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LMBR1L

**Uniprot**

Q6UX01

**Synonyms:**

489aa,55kd;

**Immunogen:**

Human LMBR1L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A