

LNX1 Antibody



PACO10275

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LNX1

Uniprot

Q8TBB1

Synonyms:

ligand of numb-protein X 1; LNX1; LNX; MPDZ; PDZRN2 ;

Immunogen:

Human LNX1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A