

LPCAT2 Antibody



PACO10280

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

LPCAT2

Uniprot

Q7L5N7

Synonyms:

lysophosphatidylcholine acyltransferase 2; LPCAT2; AYTL1; DKFZp686H22112; FLJ20481; LysoPAFAT ;

Immunogen:

Human LPCAT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A