

LRRC8D Antibody



PACO10294

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LRRC8D

Uniprot

Q7L1W4

Synonyms:

leucine rich repeat containing 8 family, member D; LRRC8D; FLJ10470; FLJ20403; LRRC5 ;

Immunogen:

Human LRRC8D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A